

**NIMBUS<sup>®</sup>**  
**PROBE**

**NIMBUS<sup>®</sup>**  
**PROBE**

**NIMBUS<sup>®</sup>**  
**PROBE**

**POST MARKET SURVEILLANCE**

Stratus Medical encourage feedback and intelligence from the clinical performance of their Electrodes to enable them to create new innovative ideas and assist with continual improvement.



Designed and manufactured by:-  
**Minta Medical Limited**  
Caddick Road,  
Knowsley Business Park,  
Knowsley, Prescot,  
Merseyside, UK,  
L34 9HP

Tel : +44 (0) 151 548 6818  
Fax : +44 (0) 151 548 5578  
[www.mintamedical.co.uk](http://www.mintamedical.co.uk)

Distributed by:  
**Stratus Medical**  
32303 FM 2978  
Magnolia, TX 77354  
USA  
[www.stratusmedical.com](http://www.stratusmedical.com)



**Radimed GmbH**  
Lothringer Strasse 36b  
44805 Bochum  
Germany  
+49 234 89 00 29-1

**DISPOSABLE RF THERMOCOUPLE  
ELECTRODE**

*For connection to:  
Neurotherm, Radionics, Baylis (Halyard),  
Cosman, Smith & Nephew, stryker,  
Owl Diros, Top and Apro Korea  
RF Lesion Generators with a  
Stratus Medical Intermediate (Adapter)  
Cable.*

Part Number

NDP-50-27-SS-CE  
NDP-100-27-SS-CE  
NDP-150-27-SS-CE

Warranty /Guarantee

This Stratus Medical disposable Thermocouple Electrode is warranted for single use only up to the use-by date shown on the pack if it has been diligently and carefully handled. It shall be observed that, due to the biological differences of the persons to be treated, no product is always absolutely effective under all environmental conditions and circumstances. Stratus Medical has no influence on the application of the product, on the diagnosis of the patient and on the handling of the product, outside of the Company. Stratus Medical can neither guarantee a beneficial or a complication-free application of the product. Stratus Medical therefore assumes no liability for damages or costs. Stratus Medical will replace products showing a deficiency which is to be represented by Stratus Medical. Employees of Stratus Medical shall not be authorised to amend the aforementioned conditions, to extend liability, or to accept or consent to additional product-related obligations.

**Instructions for use**

Caution

This product is supplied sterile for single use only. Inspect all packs before use and if the packaging is damaged do not use

Warnings

It is important that the Operating Instructions, supplied with the Lesion Generator, be read and understood before use.

[www.stratusmedical.com](http://www.stratusmedical.com)

[www.stratusmedical.com](http://www.stratusmedical.com)

#### Indications for use

Stratus Medical Gabi-L disposable RF Thermocouple Electrodes are intended for use in the lesioning of neural tissue.

#### Caution (before use)

Inspect the Thermocouple Electrode, Cable and Connector for any signs of damage.  
Do not use if any damage is observed.

#### Electrical Connections

This Thermocouple Electrode is designed for use with a RF Lesion Generator and requires the following Intermediate Cable.

DPC-NS for use with a Neurotherm and Top RF Lesion Generator.

DPC-RS for use with a Radionics RF Lesion Generator RFG-3C or RFG-3C Plus.

DPC-CS for use with a Cosman RF Lesion Generator RFG1A, RFG1B and G4.

DPC-B4 or B14 for use with a Baylis RF Lesion Generator PMG-115.

DPC-SNS for use with a Smith & Nephew RF Lesion Generator Electrothermal 20S.

DPC-ST-ALL for use with a Stryker Multigen RF Lesion Generator.

DPC-OWA for use with a Owl Diros RF Lesion Generator URF-3AP.

DPC-AK for use with a Apro Korea Coathem RF Lesion Generator for AK-F50.

After connection and before use ensure that the readings of the temperature and impedance on the Generator are within the normal range.

#### Warning

Stratus Medical RF Electrodes **MUST NOT BE BENT** when in use. Bending will seriously affect the normal characteristics and technical factors of the Electrodes.

#### Disposable Cannula

This Thermocouple Electrode has been designed to be used with a RF disposable cannula. It does not have an adjustable collar as the tolerances of manufacture make this unnecessary.  
Ensure that the correct length of disposable cannula is used.

#### Thermocouple Electrode tip position

The tip of the Thermocouple Electrode, when in the cannula, should be visible to the eye within the non-insulated tip of the cannula.  
This should always be checked with each cannula prior to use.

#### Important

Insertion of the cannula through the skin and advancement to the correct approximate site should always be done with the stylet in place in the cannula. Only the final adjustments should be done with the thermocouple in the cannula.

#### Warning

This Stratus Medical RF Thermocouple Electrode is for single use only. Do not re-use, do not re-sterilise. Re-sterilisation will seriously affect the normal characteristics and technical factors.  
The connector will not withstand multiple use.

Products are free of Latex.



Catalogue number



Serial number



Consult instructions for use



Temperature limitation



Manufacture



Date of manufacture



Batch Code / Lot Number



Sterilized using Ethylene Oxide



Use by



Do not reuse



Do not use if package is damaged

[www.stratusmedical.com](http://www.stratusmedical.com)

© Minto Medical Limited Minto Gabi-L STRATUS : Iss. 0, 05.02.2021